

TECHNICAL DATA

KLYO MEDICAL SYSTEMS, INC., Pediatric Stretcher, Model CRS10P



The KLYO MEDICAL SYSTEMS, INC. stretchers are designed to function safely with a 551 lbs (250 Kg) patient when using the pediatric stretcher CRS10P. Designed for the pediatric sector to use as transport trolley or ICU bed. Visible rails on all 4 sides ensures the best care for babies and children.

1464 NW 82nd Ave. – Doral, FL – 33126
www.klyomedical.com – sales@klyomedical.com
Tel. 1 (305) 330-5025

* Due to technical advances KLYO Medical Systems, INC., reserves the right to modify technical specifications without prior notice.

*The images shown are only for demonstrative purposes.



TECHNICAL DATA

KLYO MEDICAL SYSTEMS, INC., Pediatric Stretcher, Model CRS10P

Model	CRS10P Pediatric Stretcher
Dimension and angles have 5% tolerance.	
Specifications	
Platform Dimension (Overall Dimension)	860 x 1600 mm (970 x 1760 mm)
Adjustable Height	750~920 mm
Back Adjustment	0~50° by manual ratchet
Foot Adjustment	0~50° by manual ratchet
Trendelenburg	0~9° by foot switch
Reverse Trendelenburg	0~9° by foot switch
Hydraulic Column	Dual hydraulic pump system
Sleeping Surface	3 sections, detachable ABS material
Loading Capacity	Max. patient weight: 250 kg (551 lbs) Safety working weight: 330 kg (727.5 lbs)
Functional Equipment	
Side Rail Type	Stepless full length side rail (side rail can stop in any position during pathway). On 4 sides for 360° care. Safety switch by 2 hand operation only for lower rail.
Bed Caster	Four 8" caster; Caster control with central locking and steering system
5 th Central Wheel	One 6" caster; retractable wheel
Foldable Push Bar	4 units (front x 2, rear x 2); P shape with foam cover
IV Pole Socket	4 IV pole sockets on 4 corners
Corner Bumper	4 corners
Mattress	10 cm; waterproof
Chaise	ABS material; With oxygen holder and storage compartment

ISO 9001:2015 Certified

ISO 13485: 2016 / EN ISO 13485:2016 Certified

FDA Registered

1464 NW 82nd Ave. – Doral, FL – 33126
www.klyomedical.com – sales@klyomedical.com
 Tel. 1 (305) 330-5025

* Due to technical advances KLYO Medical Systems, INC., reserves the right to modify technical specifications without prior notice.

*The images shown are only for demonstrative purposes.