



## TECHNICAL DATA

### KLYO MEDICAL SYSTEMS, INC., Advanced Trauma Stretcher, Model CRS30



The KLYO MEDICAL SYSTEMS, INC. stretchers are designed to function safely with a 661 lbs (300 Kg) patient when using the advanced trauma stretcher CRS30. For disabled, fractured and trauma patients to transport in hospital, a built-in X-ray cassette prevents unnecessary moving and reduces risk of injury.

1464 NW 82<sup>nd</sup> Ave. – Doral, FL – 33126  
[www.klyomedical.com](http://www.klyomedical.com) – [sales@klyomedical.com](mailto:sales@klyomedical.com)  
Tel. 1 (305) 330-5025

\* Due to technical advances KLYO Medical Systems, INC., reserves the right to modify technical specifications without prior notice.

\*The images shown are only for demonstrative purposes.

DS 018-00-05/March/2021



## TECHNICAL DATA

### KLYO MEDICAL SYSTEMS, INC., Advanced Trauma Stretcher, Model CRS30

Model	CRS30 Advanced Trauma Stretcher
Dimension and angles have 5% tolerance.	
Specifications	
Dimension Platform (Dimension Overall)	619 x 1885 mm (795 x 2100 mm) – 24.37 x 74.21 in (31.30 x 82.68 in)
Adjustable Height	580~930 mm (22.84-36.61 in)
Back Adjustment	0~90° by double gas-spring
Knee Adjustment	0~40° by double gas-spring
Trendelenburg	0~18° by foot switch
Reverse Trendelenburg	0~18° by foot switch
Hydraulic Column	Dual hydraulic pump system
Sleeping Surface	4 sections, Phenolic material (X-ray translucent)
Loading Capacity	Max. patient weight: 300 kg (661 lbs) Safety working weight: 380 kg (838 lbs)
Functional Equipment	
Side Rail Type	Standard 3/4 length; Aluminum shaft with 6 units SS tube
Bed Caster	8" caster x 4; Castor control with central locking and steering system
5 <sup>th</sup> Central Wheel	One 6" caster; retractable wheel
Foldable Push Bar	4 units (front x 2, rear x 2); P shape with foam cover
Urinal Hook	Both sides
Additional Socket	4 units (front x 2, rear x 2)
IV Pole And Socket	2 IV pole units attached on diagonal corner
Side Bumper	Surrounded with plastic bumper on long sides
Corner Bumper	ABS bumper with 3" wheel protection
Mattress	7 cm; waterproof
Chaise	ABS material; With oxygen holder and storage compartment
X-Ray Cassette	Wheel operation, film cassette under platform

ISO 9001:2015 Certified

ISO 13485: 2016 / EN ISO 13485:2016 Certified

FDA Registered

1464 NW 82<sup>nd</sup> Ave. – Doral, FL – 33126  
[www.klyomedical.com](http://www.klyomedical.com) – [sales@klyomedical.com](mailto:sales@klyomedical.com)  
 Tel. 1 (305) 330-5025

\* Due to technical advances KLYO Medical Systems, INC., reserves the right to modify technical specifications without prior notice.

\*The images shown are only for demonstrative purposes.

DS 018-00-05/March/2021