




TECHNICAL DATA

KLYO MEDICAL SYSTEMS, INC., Surgical Table, Model KL10



The KLYO MEDICAL SYSTEMS, INC. tables are designed to function safely with a 792 lbs (360 Kg) patient in the mechanical model KL10. All movements are soft and silent, and the table offers a longitudinal displacement of 17.32" (44 cm), to facilitate C- Arm access. The table is radiolucent in all sections, which enables a view of the entire anatomy.

1464 NW 82nd Ave. – Doral, FL – 33126
www.klyomedical.com – sales@klyomedical.com
Tel. 1 (305) 330-5025

* Due to technical advances KLYO Medical Systems, INC., reserves the right to modify technical specifications without prior notice.

*The images shown are only for demonstrative purposes.

DS 001-01-6/29/2018



TECHNICAL DATA

KLYO MEDICAL SYSTEMS, INC., Surgical Table, Model KL10

Technical Data		KL10
Movements		Manual (mechanical)
Transfer		3" Antistatic wheels (7.6 cm)
Brake System		Manual (mechanical)
Manual Movements	Column elevation and descent	19.68" to 78.74" (80 – 100 cm)
	Removable head section	Adjustable from +45°/-45°
	Removable, independent leg plates	Continuously adjustable: UP: +38° ± 1° DOWN: -90° ± 1° Scissor movement from 0° to 180°
	Kidney elevator mechanism	4" (11 cm)
	Fowler (Back section motion)	UP: +70° ± 1° DOWN: -35° ± 1°
	Trendelenburg	25° ± 1°
	Reverse Trendelenburg	25° ± 1°
	Left Tilt	15° ± 1°
	Right Tilt	15° ± 1°
	Longitudinal Displacement	17.32" (44 cm)
	Flex-Reflex	yes
Patient Orientation Selection		Yes, through manual Exchange of back/head section with leg section.
Maximum patient weight		792 lbs (360 Kg)
Base		Stainless Steel
Base cover material		ABS or stainless steel
Column		Stainless Steel
Surface Structure		Stainless Steel
Side Rails		Stainless Steel
Table Top	Material	Fenolite, completely radiolucent
	Number of Sections	5 to 7 ▪ Head ▪ Separate leg plates ▪ Pelvis ▪ Kidney elevator ▪ Complete back (Optional: bipartite)
Length		78.74" (200 cm)
Width (including side rails)		19.68" (50 cm)
Pads	Material	Viscoelastic, removable, independent, antistatic, seamless.
	Thickness	2" (5.08 cm) optional 3" (7.62 cm)
	Radiolucent	yes

1464 NW 82nd Ave. – Doral, FL – 33126
www.klyomedical.com – sales@klyomedical.com
Tel. 1 (305) 330-5025

* Due to technical advances KLYO Medical Systems, INC., reserves the right to modify technical specifications without prior notice.

*The images shown are only for demonstrative purposes.

DS 001-01-6/29/2018



TECHNICAL DATA

KLYO MEDICAL SYSTEMS, INC., Surgical Table, Model KL10



QUALITY CERTIFICATIONS

Amtivo (USA) Inc.

SO 9001:2015 Certified ISO 13485: 2016

EN ISO 13485:2016 Certified

FDA Certified

1464 NW 82nd Ave. – Doral, FL – 33126
www.klyomedical.com – sales@klyomedical.com
Tel. 1 (305) 330-5025

* Due to technical advances KLYO Medical Systems, INC., reserves the right to modify technical specifications without prior notice.

*The images shown are only for demonstrative purposes.

DS 001-01-6/29/2018